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Tampon super-absorbant

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Indications :

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Dispositif Médical

Device Identification

Super-absorbent pad
Presentation: by 25

Indications:

Gels body fluids in a collection bin, waste container or various bedpan-type containers. It can also be used to increase the absorption capacity of CareBag® bags, if required. This device is particularly intended for healthcare professionals.

Composition:

Cellulose based support + sodium polyacrylate.

Performance & mechanism of action:

Upon contact with the super-absorbent SAP component (contained in the pad) the excreted matter gel. SAPs are polymers that absorb and retain large volumes of liquid compared to their own mass. They remain insoluble in the liquid they absorb (thanks to their cross-linked structure). The absorption capacity of each pad is 450ml of 0.9% saline solution.

Instructions for use:

1. Open the pack and take a super-absorbent pad,
2. Place the super-absorbent pad in the desired area,
3. Wait for liquids to turn to gel,
4. Dispose of in appropriate place.

When using with CAREBAG® products, it is recommended to add pads according to the volume of additional fluids to be collected:

CareBag® range	Bedpan Liner / Commode Liner / Toilet Liner	Male Urinal / Vomit bag
Maximum volume absorption	3 L	1,5 L
Number of Gelmax® Pads	6	3 - 4

Warnings, precautions for use and contraindications:

- For single use. In case of reuse, there is a risk of insufficient or non-existent absorption, and therefore of potential contamination and discomfort for the patient,
- Do not use on damaged skin;
- Do not use this product to gel the following chemical products: strong acids (concentrations > 3%), strong bases (concentrations > 10%) and strong oxidizing agents,
- Keep away from children,
- In case the pad comes into contact with the eye(s) or mucous membranes, rinse with plenty of tap water. In case of persistent discomfort, contact your doctor,
- In case of ingestion, contact a doctor immediately.

Side effects:

Any serious incident with the medical device should be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Storage & Disposal:

Product must be kept dry and shall be stored away from sunlight. Dispose of in household waste, or according to the waste management protocol in place at the establishment.

Publication date : 01/2024

Single use

Keep dry

Keep away from the sunlight

Consult the instructions for use

Medical Device